

Product brochure

07.08.2026

KERN 304-07

Weight set, E1, 1 g - 2 kg, Knob, stainless steel polished, in wooden case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-070-100	Wood weight case 313-070-100 (for nominal values 1 mg - 2 kg for classes E1+E2+F1 for design Cylindrical/compact)

SERVICES

Model	Description
920-715	Marking (Last digits of serial number) KERN 920-715
921-220	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-220
962-220R	DAkKS Recalibration 962-220R
963-220	DAkKS calibration certificate 963-220

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days