

Product brochure

07.08.2026

# KERN 307-09

Individual weight, OIML Class E1, 500 g, Knob, stainless steel polished

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com

## PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

## PICTOGRAMS

Standard

Option



## ACCESSORIES

| Model       | Description                                                                                                         |
|-------------|---------------------------------------------------------------------------------------------------------------------|
| 315-245     | Tweezers 315-245 (stainless steel   250 mm)                                                                         |
| 317-090-100 | Wood weight case 317-090-100 (for nominal values 500 g   for classes E1+E2+F1   for design Cylindrical/compact)     |
| 317-090-400 | Plastic weight case 317-090-400 (for nominal values 500 g   for classes E1, E2   for design Cylindrical/compact)    |
| 317-090-600 | Aluminium weight case 317-090-600 (for nominal values 500 g   for classes E1 - M3   for design Cylindrical/compact) |
| 317-280     | Gloves 317-280 (cotton, 1 pair)                                                                                     |

## SERVICES

| Model    | Description                                                         |
|----------|---------------------------------------------------------------------|
| 920-350  | Marking by laser engraving (individual marking) KERN 920-350        |
| 920-351  | Marking of test weights (last digits of serial number) KERN 920-351 |
| 921-239  | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-239             |
| 962-239R | DAkKS Recalibration 962-239R                                        |
| 963-239  | DAkKS calibration certificate 963-239                               |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days