

Product brochure

07.08.2026

KERN 307-11

Individual weight, OIML Class E1, 1 kg, Knob, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-245	Tweezers 315-245 (stainless steel 250 mm)
317-110-100	Wood weight case 317-110-100 (for nominal values 1 kg for classes E1+E2+F1 for design Cylindrical/compact)
317-110-400	Plastic weight case 317-110-400 (for nominal values 1 kg for classes E2 for design Cylindrical/compact)
317-110-600	Aluminium weight case 317-110-600 (for nominal values 1 kg for classes E1 - M3 for design Cylindrical/compact)
317-280	Gloves 317-280 (cotton, 1 pair)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-241	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-241
962-241R	DAkkS Recalibration 962-241R
963-241	DAkkS calibration certificate 963-241

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days