

Product brochure

07.08.2026

KERN 307-16

Individual weight, OIML Class E1, 50 kg, Knob, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
317-160-100	Wood weight case 317-160-100 (for nominal values 50 kg for classes E1+E2+F1 for design Cylindrical/compact)
317-160-600	Aluminium weight case 317-160-600 (for nominal values 50 kg for classes E1 - M3 for design Cylindrical/compact)
317-290	Gloves 317-290 (leather, cotton, 1 pair)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-246	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-246
962-246R	DAkKS Recalibration 962-246R
963-246	DAkKS calibration certificate 963-246

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days