

# KERN 307-16

Individual weight, OIML Class E1, 50 kg, Knob, Stainless steel polished



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	307

## Measuring System

Conformity assessment possible	✓
Test weight - volume (cm <sup>3</sup> )	6250
Test weight - volume uncertainty (cm <sup>3</sup> )	78,125
Test weight - Density (kg/m <sup>3</sup> )	8000
Test weight - Density uncertainty (kg/m <sup>3</sup> )	100
Test weight - Density [Min] (OIML)	7934
Test weight - Density [Max] (OIML)	8067
Test weight - nominal value	50 kg
OIML-compliant	✓
Weight OIML class	E1
Permissible tolerance according to OIML (+/-)	25 mg

## Construction

Shape weight	Knob
Dimensions (Ø×H)	178×298 mm
Material	Stainless steel polished

## Services (optional)

Article number for DAkkS calibration	963-246
Article number for guaranteed plus-deviation	921-246

## Packing & Shipping

Delivery time	2 d
Dimensions packaging (W×D×H)	178×178×298 mm
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg

## Optional cases

Case - Aluminium article number	317-160-600
Case - Wood article number	317-160-100

## Pictograms

### STANDARD



### OPTION

