

# KERN 308-36

**KERN**

Milligram weights, E1, 50 mg, Wire, stainless steel polished



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	308

## Measuring System

Conformity assessment possible	✓
Test weight - volume (cm <sup>3</sup> )	0,0063
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0001
Test weight - Density (kg/m <sup>3</sup> )	8000
Test weight - Density uncertainty (kg/m <sup>3</sup> )	160
Test weight - Density [Min] (OIML)	3400
Test weight - Density [Max] (OIML)	100000
Test weight - nominal value	50 mg
OIML-compliant	✓
Weight OIML class	E1
Permissible tolerance according to OIML (+/-)	0,004 mg

## Construction

Shape weight	Wire
Dimensions weight (W×D×H)	11×11×3 mm
Material	stainless steel polished

## Services (optional)

Article number for DAkkS calibration	962-256
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	11×11×3 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

## Pictograms

### STANDARD

