

KERN 308-37

KERN

Milligram weights, E1, 100 mg, Wire, stainless steel polished



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	308

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	0,0125
Test weight - volume uncertainty (cm ³)	0,0003
Test weight - Density (kg/m ³)	8000
Test weight - Density uncertainty (kg/m ³)	160
Test weight - Density [Min] (OIML)	4400
Test weight - Density [Max] (OIML)	100000
Test weight - nominal value	100 mg
OIML-compliant	✓
Weight OIML class	E1
Permissible tolerance according to OIML (+/-)	0,005 mg

Construction

Shape weight	Wire
Dimensions weight (W×D×H)	9×9×3 mm
Material	stainless steel polished

Services (optional)

Article number for DAkkS calibration	962-257
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	9×9×3 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

Pictograms

STANDARD

