

Product brochure

08.08.2026

KERN 308-426

Set of milligram weights, E1, 1 mg - 500 mg, Wire, stainless steel polished, in aluminium case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-247	Tweezers 315-247 (stainless steel, plastic 136 mm)

SERVICES

Model	Description
920-701	Marking (Last digits of serial number) KERN 920-701
962-250	DAkKS calibration certificate 962-250
962-250R	DAkKS Recalibration 962-250R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days