

KERN 313-010-600

KERN

Aluminium protective box for weight sets (Fracture gram, 1 mg - 500 mg)



Pictograms

STANDARD



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Aluminium protective box

Measuring System

Case suitable for OIML classes	E1 - M2
Case suitable for nominal weight values	1 mg - 500 mg
Suitable weight shape for case	Fracture gram

Construction

Dimension (W×D×H)	165×158×37 mm
Material	aluminium
Case - Smallest possible weight	1 mg
Case - Highest possible weight	500 mg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×158×37 mm
Net weight	0,3 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,30 kg
Shipping weight	0,3 kg