

# KERN 313-030-100

Wooden box for weight sets (Cylindrical/compact, 1 mg - 100 g)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	E1+E2
Case suitable for nominal weight values	1 mg - 100 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	152×152×66 mm
Material	beech
Case - Smallest possible weight	1 mg
Case - Highest possible weight	100 g

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	152×152×66 mm
Net weight	0,68 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	0,70 kg
Shipping weight	0,68 kg

## Pictograms

### STANDARD

