

# KERN 313-040-600

**KERN**

Aluminium protective box for weight sets (Cylindrical/compact, 1 mg - 200 g)



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Aluminium protective box

## Measuring System

Case suitable for OIML classes	E1 - M1
Case suitable for nominal weight values	1 mg - 200 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	190×130×80 mm
Material	aluminium
Case - Smallest possible weight	1 mg
Case - Highest possible weight	200 g

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	190×130×80 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg