

Product brochure

07.08.2026

# KERN 313-044

Weight set, E2, 1 mg - 200 g, Knob, stainless steel polished, in plastic case


ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

| Model       | Description                                                                                                             |
|-------------|-------------------------------------------------------------------------------------------------------------------------|
| 313-052-400 | Plastic carrying case 313-052-400 (for nominal values bis 500 g   for classes E1 - M3   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                          |
|----------|------------------------------------------------------|
| 920-704  | Marking (Last digits of serial number) KERN 920-704  |
| 921-303  | Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-303 |
| 962-303  | DAkKS calibration certificate 962-303                |
| 962-303R | DAkKS Recalibration 962-303R                         |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days