

KERN 313-050-100

KERN

Wooden box for weight sets (Cylindrical/compact, 1 mg - 500 g)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	1 mg - 500 g
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	230×180×110 mm
Material	beech
Case - Smallest possible weight	1 mg
Case - Highest possible weight	500 g

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	230×180×110 mm
Net weight	2,1 kg
Shipping method	Parcel service
Net weight approx.	2,2 kg
Gross weight approx.	2,6 kg
Shipping weight	2,5 kg

Pictograms

STANDARD

