

KERN 313-050-400

KERN

Plastic box for weight sets (Cylindrical/compact, bis 500 g)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Plastic box

Measuring System

Case suitable for OIML classes	E1 - M3
Case suitable for nominal weight values	bis 500 g
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	250×210×75 mm
Material	plastic
Case - Smallest possible weight	1 mg
Case - Highest possible weight	500 g

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	0,35 kg
Shipping method	Parcel service
Net weight approx.	0,35 kg
Gross weight approx.	0,35 kg
Shipping weight	0,8 kg

Pictograms

STANDARD

