

KERN 313-080-100

Wooden box for weight sets (Cylindrical/compact, 1 mg - 5 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	1 mg - 5 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	420×200×200 mm
Material	beech
Case - Smallest possible weight	1 mg
Case - Highest possible weight	5 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	420×200×200 mm
Net weight	6,7 kg
Shipping method	Parcel service
Net weight approx.	7 kg
Gross weight approx.	8 kg
Shipping weight	7,3 kg

Pictograms

STANDARD

