

Product brochure

08.08.2026

KERN 313-086

Weight set, E2, 1 mg - 5 kg, Knob, stainless steel polished, in aluminium case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-080-600	Aluminium weight case 313-080-600 (for nominal values 1 mg - 5 kg for classes E1 - M1 for design Cylindrical/compact)

SERVICES

Model	Description
920-708	Marking (Last digits of serial number) KERN 920-708
921-307	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-307
962-307	DAkKS calibration certificate 962-307
962-307R	DAkKS Recalibration 962-307R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days