

Product brochure

07.08.2026

KERN 313-09

Weight set, E2, 1 mg - 10 kg, Knob, stainless steel polished, in wooden case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



FEATURES

- The milligram weights contained in this weight set (denomination 1 mg - 500 mg) are in a separate plastic case

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
313-090-100	Wood weight case 313-090-100 (for nominal values 1 mg - 10 kg for classes E1+E2+F1 for design Cylindrical/compact)

SERVICES

Model	Description
920-709	Marking (Last digits of serial number) KERN 920-709
921-308	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-308
962-308	DAkKS calibration certificate 962-308
962-308R	DAkKS Recalibration 962-308R

KERN Pictogram Overview



Pallet shipment: The duration of internal product provision are two days



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days