

Product brochure

07.08.2026

KERN 313-090-100

Wooden box for weight sets (Cylindrical/compact, 1 mg - 10 kg)

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
317-280	Gloves 317-280 (cotton, 1 pair)
318-270	Brush 318-270
328-222-410	Weight case 328-222-410 (for nominal values 1 mg - 5 g for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day