

# KERN 313-090-100

Wooden box for weight sets (Cylindrical/compact, 1 mg - 10 kg)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	1 mg - 10 kg
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	540×200×240 mm
Material	beech
Case - Smallest possible weight	1 mg
Case - Highest possible weight	10 kg

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	540×200×240 mm
Net weight	9,5 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	11 kg
Shipping weight	10,5 kg

## Pictograms

### STANDARD

