

Product brochure

07.08.2026

KERN 314-024

Weight set, E2, 1 g - 50 g, Knob, stainless steel polished, in plastic case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-052-400	Plastic carrying case 313-052-400 (for nominal values bis 500 g for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-710	Marking (Last digits of serial number) KERN 920-710
921-315	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-315
962-315	DAkkS calibration certificate 962-315
962-315R	DAkkS Recalibration 962-315R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days