

# KERN 314-026

# KERN

Weight set, E2, 1 g - 50 g, Knob, stainless steel polished, in aluminium case



Shipping weight 0,8 kg

## Optional cases

Case - Aluminium article number 314-020-600

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	314

## Measuring System

Conformity assessment possible	✓
Test weight - nominal value	1 g - 50 g
OIML-compliant	✓
Weight OIML class	E2

## Construction

Dimension (W×D×H)	165×120×75 mm
Shape weight	Knob
Material	stainless steel polished
Case of weight set	in aluminium case

## Services (optional)

Article number for conformity assessment	952-315
Article number for DAkkS calibration	962-315
Article number for guaranteed plus-deviation	921-315

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×120×75 mm
Shipping method	Parcel service
Net weight approx.	0,55 kg
Gross weight approx.	0,80 kg