

Product brochure

07.08.2026

# KERN 314-044

Weight set, E2, 1 g - 200 g, Knob, stainless steel polished, in plastic case

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

☎ +49 7433 9933-0  
✉ [info@kern-sohn.com](mailto:info@kern-sohn.com)  
🌐 [www.kern-sohn.com](http://www.kern-sohn.com)



PICTOGRAMS

Standard

Option

ACCESSORIES

| Model       | Description                                                                                                             |
|-------------|-------------------------------------------------------------------------------------------------------------------------|
| 313-052-400 | Plastic carrying case 313-052-400 (for nominal values bis 500 g   for classes E1 - M3   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                             |
|----------|---------------------------------------------------------|
| 920-712  | Marking (Last digits of serial number) KERN 920-712     |
| 921-317  | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-317 |
| 962-317  | DAkKS calibration certificate 962-317                   |
| 962-317R | DAkKS Recalibration 962-317R                            |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days