

Product brochure

08.08.2026

KERN 314-046

Weight set, E2, 1 g - 200 g, Knob, stainless steel polished, in aluminium case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-040-600	Aluminium weight case 313-040-600 (for nominal values 1 mg - 200 g for classes E1 - M1 for design Cylindrical/compact)

SERVICES

Model	Description
920-712	Marking (Last digits of serial number) KERN 920-712
921-317	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-317
962-317	DAKKS calibration certificate 962-317
962-317R	DAKKS Recalibration 962-317R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days