

Product brochure

07.08.2026

# KERN 314-05

Weight set, E2, 1 g - 500 g, Knob, stainless steel polished, in wooden case


ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

| Model       | Description                                                                                                            |
|-------------|------------------------------------------------------------------------------------------------------------------------|
| 313-050-100 | Wood weight case 313-050-100 (for nominal values 1 mg - 500 g   for classes E1+E2+F1   for design Cylindrical/compact) |

SERVICES

| Model    | Description                                             |
|----------|---------------------------------------------------------|
| 920-713  | Marking (Last digits of serial number) KERN 920-713     |
| 921-318  | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-318 |
| 962-318  | DAkkS calibration certificate 962-318                   |
| 962-318R | DAkkS Recalibration 962-318R                            |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days