

KERN 314-056

KERN

Weight set, E2, 1 g - 500 g, Knob, stainless steel polished, in aluminium case



Shipping weight 2 kg

Optional cases

Case - Aluminium article number 313-050-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	314

Measuring System

Conformity assessment possible	✓
Test weight - nominal value	1 g - 500 g
OIML-compliant	✓
Weight OIML class	E2

Construction

Dimension (W×D×H)	210×140×98 mm
Shape weight	Knob
Material	stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for conformity assessment	952-318
Article number for DAkkS calibration	962-318
Article number for guaranteed plus-deviation	921-318

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	210×140×98 mm
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	2,0 kg