

KERN 314-06

KERN

Weight set, E2, 1 g - 1 kg, Knob, stainless steel polished, in wooden case



Shipping weight 4,21 kg

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	314

Measuring System

Test weight - nominal value	1 g - 1 kg
OIML-compliant	✓
Weight OIML class	E2

Construction

Dimension (W×D×H)	230×180×110 mm
Shape weight	Knob
Material	stainless steel polished
Number of weights	13
Case of weight set	in wooden case

Services (optional)

Article number for conformity assessment	952-319
Article number for DAkkS calibration	962-319
Article number for guaranteed plus-deviation	921-319

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	230×180×110 mm
Shipping method	Parcel service
Net weight approx.	4,4 kg
Gross weight approx.	4,4 kg