

Product brochure

08.08.2026

KERN 314-066

Weight set, E2, 1 g - 1 kg, Knob, stainless steel polished, in aluminium case

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-060-600	Aluminium weight case 313-060-600 (for nominal values 1 mg - 1 kg for classes E1 - M1 for design Cylindrical/compact)

SERVICES

Model	Description
920-714	Marking (Last digits of serial number) KERN 920-714
921-319	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-319
962-319	DAkKS calibration certificate 962-319
962-319R	DAkKS Recalibration 962-319R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days