

Product brochure

08.08.2026

KERN 314-076

Weight set, E2, 1 g - 2 kg, Knob, stainless steel polished, in aluminium case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



PICTOGRAMS

Standard

Option

ACCESSORIES

Model	Description
313-070-600	Aluminium weight case 313-070-600 (for nominal values 1 mg - 2 kg for classes E1 - M1 for design Cylindrical/compact)

SERVICES

Model	Description
920-715	Marking (Last digits of serial number) KERN 920-715
921-320	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-320
962-320	DAkkS calibration certificate 962-320
962-320R	DAkkS Recalibration 962-320R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days