

# KERN 314-076

**KERN**

Weight set, E2, 1 g - 2 kg, Knob, stainless steel polished, in aluminium case



Shipping weight 5,55 kg

## Optional cases

Case - Aluminium article number 313-070-600

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E2
Product family	314

## Measuring System

Conformity assessment possible	✓
Test weight - nominal value	1 g - 2 kg
OIML-compliant	✓
Weight OIML class	E2

## Construction

Dimension (W×D×H)	278×180×138 mm
Shape weight	Knob
Material	stainless steel polished
Case of weight set	in aluminium case

## Services (optional)

Article number for conformity assessment	952-320
Article number for DAkkS calibration	962-320
Article number for guaranteed plus-deviation	921-320

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	278×180×138 mm
Shipping method	Parcel service
Net weight approx.	6 kg
Gross weight approx.	6 kg