

Product brochure

07.08.2026

KERN 314-084

Weight set, E2, 1 g - 5 kg, Knob, stainless steel polished, in plastic case


ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



PICTOGRAMS

Standard

Option





ACCESSORIES

Model	Description
313-082-400	Plastic carrying case 313-082-400 (for nominal values bis 5 kg for classes E1 - M3 for design Cylindrical/compact)

SERVICES

Model	Description
920-716	Marking (Last digits of serial number) KERN 920-716
921-321	Guaranteed Plus-deviation (> +1/3 OIML) KERN 921-321
962-321	DAkkS calibration certificate 962-321
962-321R	DAkkS Recalibration 962-321R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days