

# KERN 315-070-100

Wooden box for individual weight sets



## Category

|                  |               |
|------------------|---------------|
| Brand            | KERN          |
| Product category | Cases & Boxes |
| Product group    | Wooden box    |

## Measuring System

|   |                     |
|---|---------------------|
| Case suitable for OIML classes          | E1+E2+F1            |
| Case suitable for nominal weight values | universal bis 2 kg  |
| Suitable weight shape for case          | Cylindrical/compact |

## Construction

|                                 |                |
|---------------------------------|----------------|
| Dimension (W×D×H)               | 305×180×140 mm |
| Material                        | beech          |
| Case - Smallest possible weight | 1 mg           |
| Case - Highest possible weight  | 2 kg           |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Delivery time                | 1 d            |
| Dimensions packaging (W×D×H) | 310×180×140 mm |
| Net weight                   | 2,9 kg         |
| Shipping method              | Parcel service |
| Net weight approx.           | 3,0 kg         |
| Gross weight approx.         | 3,0 kg         |
| Shipping weight              | 2,9 kg         |

## Pictograms

### STANDARD

