

Product brochure

07.08.2026

KERN 315-080-100

Wooden box for individual weight sets

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

FEATURES

- Configure your set of weights entirely according to your requirements
- KERN manufactures a tailor-made wooden box or alu or plastic case for you – shaped and upholstered according to your requirements
- Please note: These universal cases are customised and therefore not eligible for return

PLEASE NOTE FOR YOUR ORDER

- Please note: These universal cases are customised and therefore not eligible for return

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
315-242	Tweezers 315-242 (polyamide 100 mm)
315-243	Tweezers 315-243 (stainless steel 105 mm)
317-280	Gloves 317-280 (cotton, 1 pair)
318-270	Brush 318-270
328-222-410	Weight case 328-222-410 (for nominal values 1 mg - 5 g for classes E1 - M2 for design Cylindrical/compact)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day