


KERN 315-247



Tweezers to be able to safely grip small test weights



Category

Brand	KERN
Product category	Test weights accessories
Product group	Tweezer

Measuring System

Case suitable for OIML classes	E1 - M3
--------------------------------	---------

Construction

Material	stainless steel, plastic
----------	--------------------------

Environmental conditions

Ambient temperature [Max]	140 °C
---------------------------	--------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	135×15×10 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Pictograms

STANDARD

