

Product brochure

08.08.2026

KERN 315-308-204

Safety Set 200 g (F2); 10 g (F1)

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com



FEATURES

- These “KERN Safety Sets” each contain one test weight for testing the sensitivity, i.e. the correct adjustment of your balance as well as a small test weight to test the lower end of the weighing range, the so-called minimum sample weight
- As an option, the “KERN Safety Set” has space for another test weight, for testing your balance at a weight which is relevant for you
- Useful accessories which have been selected to suit that particular “KERN Safety Set”, such as, for example, special gloves, tweezers, weight grips, brushes, etc. will assist you in handling your test weights properly
- Stored in the practical protective case next to your balance, you can check and ensure the high precision of your balance at any time

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
921-438	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-438
962-400-600	DAkkS calibration certificate 962-400-600
962-600-600	DAkkS calibration certificate 962-600-600

SERVICES

Model	Description
921-438	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-438
962-400-400	DAkkS calibration certificate 962-400-400

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day