

Product brochure

08.08.2026

# KERN 315-312-308

Safety Set 200 g (F2); 2000 g (F2)

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com



## FEATURES

- These “KERN Safety Sets” each contain one test weight for testing the sensitivity, i.e. the correct adjustment of your balance as well as a small test weight to test the lower end of the weighing range, the so-called minimum sample weight
- As an option, the “KERN Safety Set” has space for another test weight, for testing your balance at a weight which is relevant for you
- Useful accessories which have been selected to suit that particular “KERN Safety Set”, such as, for example, special gloves, tweezers, weight grips, brushes, etc. will assist you in handling your test weights properly
- Stored in the practical protective case next to your balance, you can check and ensure the high precision of your balance at any time

## PICTOGRAMS

Standard



## ACCESSORIES

| Model       | Description                                             |
|-------------|---------------------------------------------------------|
| 921-442     | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-442 |
| 962-400-600 | DAkkS calibration certificate 962-400-600               |
| 962-600-600 | DAkkS calibration certificate 962-600-600               |

## SERVICES

| Model       | Description                                             |
|-------------|---------------------------------------------------------|
| 921-442     | Guaranteed Plus-deviation ( $> +1/3$ OIML) KERN 921-442 |
| 962-400-400 | DAkkS calibration certificate 962-400-400               |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day