

KERN 316-05

Individual weight, OIML Class E2, 20 g, Compact form, stainless steel polished



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	316

Measuring System

Conformity assessment possible	no
Test weight - volume (cm ³)	2,5
Test weight - volume uncertainty (cm ³)	0,0313
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	7500
Test weight - Density [Max] (OIML)	8570
Test weight - nominal value	20 g
OIML-compliant	✓
Weight OIML class	E2
Permissible tolerance according to OIML (+/-)	0,08 mg

Construction

Shape weight	Compact form
Dimensions (Ø×H)	13×20 mm
Material	stainless steel polished

Services (optional)

Article number for DAkkS calibration	962-335
Article number for guaranteed plus-deviation	921-335

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	13×13×20 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	317-050-400
Case - Aluminium article number	317-050-600
Case - Wood article number	317-050-100

Pictograms

STANDARD

