

Product brochure

07.08.2026

KERN 316-09

Individual weight, OIML Class E2, 500 g, Compact form, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
315-245	Tweezers 315-245 (stainless steel 250 mm)
317-090-100	Wood weight case 317-090-100 (for nominal values 500 g for classes E1+E2+F1 for design Cylindrical/compact)
317-090-400	Plastic weight case 317-090-400 (for nominal values 500 g for classes E1, E2 for design Cylindrical/compact)
317-090-600	Aluminium weight case 317-090-600 (for nominal values 500 g for classes E1 - M3 for design Cylindrical/compact)
317-280	Gloves 317-280 (cotton, 1 pair)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-339	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-339
962-339	DAkKS calibration certificate 962-339
962-339R	DAkKS Recalibration 962-339R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days