

KERN 317-009-600

KERN

Aluminium protective box for individual weights (Fracture gram, 1 mg - 500 mg (einzel))



Pictograms

STANDARD



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Aluminium protective box

Measuring System

Case suitable for OIML classes	E1 - M3
Case suitable for nominal weight values	1 mg - 500 mg (einzel)
Suitable weight shape for case	Fracture gram

Construction

Dimension (W×D×H)	70×77×73 mm
Material	aluminium
Case - Smallest possible weight	1 mg
Case - Highest possible weight	500 mg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	77×70×73 mm
Net weight	0,3 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,30 kg
Shipping weight	0,3 kg