

Product brochure

07.08.2026

KERN 317-030-100

Wooden box for individual weights (Cylindrical/compact, 5 g)

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com

FEATURES

- Suitable for individual weights KERN 307, 316, 317, 326, 327, 329

PICTOGRAMS

Standard



SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day