

Product brochure

07.08.2026

# KERN 317-070-400

Plastic box for individual weights (Cylindrical/compact, 100 g)

ENGLISH



## CONTACT

KERN & SOHN GmbH  
Ziegelei 1  
72336 Balingen  
Germany

+49 7433 9933-0  
info@kern-sohn.com  
www.kern-sohn.com



FEATURES

- Box for individual weights

PICTOGRAMS

Standard



SERVICES

| Model   | Description                                                  |
|---------|--------------------------------------------------------------|
| 920-350 | Marking by laser engraving (individual marking) KERN 920-350 |

## KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day