

# KERN 317-090-100



Wooden box for individual weights (Cylindrical/compact, 500 g)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	500 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	70×70×90 mm
Material	beech
Case - Smallest possible weight	500 g
Case - Highest possible weight	500 g

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	70×70×90 mm
Net weight	0,34 kg
Shipping method	Parcel service
Net weight approx.	0,35 kg
Gross weight approx.	0,35 kg
Shipping weight	0,34 kg

## Pictograms

### STANDARD

