

KERN 317-120-400

KERN

Plastic box for individual weights (Cylindrical/compact, 2 kg)



Category

| | |
|------------------|---------------|
| Brand | KERN |
| Product category | Cases & Boxes |
| Product group | Plastic box |

Measuring System

| | |
|---|---------------------|
| Case suitable for OIML classes | E2 |
| Case suitable for nominal weight values | 2 kg |
| Suitable weight shape for case | Cylindrical/compact |

Construction

| | |
|---------------------------------|------------|
| Dimension (Ø×H) | 100×165 mm |
| Material | plastic |
| Case - Smallest possible weight | 2 kg |
| Case - Highest possible weight | 2 kg |

Packing & Shipping

| | |
|------------------------------|----------------|
| Delivery time | 1 d |
| Dimensions packaging (W×D×H) | 100×100×165 mm |
| Net weight | 0,02 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,05 kg |
| Gross weight approx. | 0,05 kg |
| Shipping weight | 0,3 kg |

Pictograms

STANDARD

