

# KERN 317-13



Individual weight, OIML Class E2, 5 kg, Knob, stainless steel polished



Article number for guaranteed plus-deviation	921-343
--	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	80×80×140 mm
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg

## Optional cases

Case - Plastic article number	317-130-400
Case - Aluminium article number	317-130-600
Case - Wood article number	317-130-100

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

## Measuring System

Conformity assessment possible	✓
Test weight - volume (cm <sup>3</sup> )	625
Test weight - volume uncertainty (cm <sup>3</sup> )	7,8125
Test weight - Density (kg/m <sup>3</sup> )	7950
Test weight - Density uncertainty (kg/m <sup>3</sup> )	140
Test weight - Density [Min] (OIML)	7810
Test weight - Density [Max] (OIML)	8210
Test weight - nominal value	5 kg
OIML-compliant	✓
Weight OIML class	E2
Permissible tolerance according to OIML (+/-)	8 mg

## Construction

Shape weight	Knob
Dimensions (Ø×H)	80×140 mm
Material	stainless steel polished

## Services (optional)

Article number for conformity assessment	952-343
Article number for DAkkS calibration	962-343