

Product brochure

07.08.2026

KERN 317-130-400

Plastic box for individual weights (Cylindrical/compact, 5 kg)



ENGLISH

CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

+49 7433 9933-0
info@kern-sohn.com
www.kern-sohn.com



FEATURES

- Box for individual weights

PICTOGRAMS

Standard



SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day