

Product brochure

07.08.2026

KERN 317-14

Individual weight, OIML Class E2, 10 kg, Knob, stainless steel polished

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

PLEASE NOTE FOR YOUR ORDER

- For this weight, we strongly recommend a box due to the very accurate weight class. You will find a selection in the accessories section

PICTOGRAMS

Standard

Option



ACCESSORIES

Model	Description
317-140-100	Wood weight case 317-140-100 (for nominal values 10 kg for classes E1+E2+F1 for design Cylindrical/compact)
317-140-400	Plastic weight case 317-140-400 (for nominal values 10 kg for classes E2 for design Cylindrical/compact)
317-140-600	Aluminium weight case 317-140-600 (for nominal values 10 kg for classes E1 - M3 for design Cylindrical/compact)
317-280	Gloves 317-280 (cotton, 1 pair)
317-290	Gloves 317-290 (leather, cotton, 1 pair)

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350
920-351	Marking of test weights (last digits of serial number) KERN 920-351
921-344	Guaranteed Plus-deviation ($> +1/3$ OIML) KERN 921-344
962-344	DAkKS calibration certificate 962-344
962-344R	DAkKS Recalibration 962-344R

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day



Accredited calibration (DKD): The duration of the accredited calibration is 3 working days