

KERN 317-150-100

KERN

Wooden box for individual weights (Cylindrical/compact, 20 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	20 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	200×200×315 mm
Material	beech
Case - Smallest possible weight	20 kg
Case - Highest possible weight	20 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	200×200×315 mm
Net weight	4,6 kg
Shipping method	Parcel service
Net weight approx.	4,6 kg
Gross weight approx.	4,6 kg
Shipping weight	4,6 kg

Pictograms

STANDARD

