

Product brochure

07.08.2026

KERN 317-160-100

Wooden box for individual weights (Cylindrical/compact, 50 kg)

ENGLISH



CONTACT

KERN & SOHN GmbH
Ziegelei 1
72336 Balingen
Germany

☎ +49 7433 9933-0
✉ info@kern-sohn.com
🌐 www.kern-sohn.com

FEATURES

- Suitable for individual weights KERN 307, 316, 317, 326, 327, 329

PICTOGRAMS

Standard



ACCESSORIES

Model	Description
970-104	Warranty extension (+2 years) 970-104

SERVICES

Model	Description
920-350	Marking by laser engraving (individual marking) KERN 920-350

ENGLISH

KERN Pictogram Overview



Package shipment by courier service: The duration of internal product provision is one day