

KERN 317-160-100

KERN

Wooden box for individual weights (Cylindrical/compact, 50 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	50 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	185×185×300 mm
Material	beech
Case - Smallest possible weight	50 kg
Case - Highest possible weight	50 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	185×185×300 mm
Net weight	8,5 kg
Shipping method	Parcel service
Net weight approx.	9 kg
Gross weight approx.	9 kg
Shipping weight	8,5 kg

Pictograms

STANDARD

