

KERN 318-22

KERN

Set of milligram weights, E2, 1 mg - 500 mg, Platelet, stainless steel polished, in plastic case



Shipping weight 0,1 kg

Optional cases

Case - Plastic article number 328-222-420

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	318

Measuring System

Test weight - nominal value	1 mg - 500 mg
OIML-compliant	✓
Weight OIML class	E2

Construction

Dimension (W×D×H)	136×53×22 mm
Shape weight	Platelet
Material	stainless steel polished
Number of weights	13
Case of weight set	in plastic case

Services (optional)

Article number for conformity assessment	952-350
Article number for DAkkS calibration	962-350
Article number for guaranteed plus-deviation	921-350

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	130×60×22 mm
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg