

KERN 318-35

KERN

Milligram weights, E2, 20 mg, Wire, stainless steel polished



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	318

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	0,0025
Test weight - volume uncertainty (cm ³)	0,0001
Test weight - Density (kg/m ³)	8000
Test weight - Density uncertainty (kg/m ³)	160
Test weight - nominal value	20 mg
OIML-compliant	✓
Weight OIML class	E2
Permissible tolerance according to OIML (+/-)	0,01 mg

Construction

Shape weight	Wire
Dimensions weight (W×D×H)	7×7×3 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-355
Article number for DAkkS calibration	962-355
Article number for guaranteed plus-deviation	921-355

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	7×7×3 mm
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

Pictograms

STANDARD

